

NAF  **NaturalintTMX**

FIRST AID

Dressing

Highly absorbent quality
cotton wool, encased in
a low-adherent, tubular
non-woven cover

500g



**Cushioned padding
Insulation for
Minor wound dressing**

DIRECTIONS FOR USE: *Clean hands thoroughly before and after application. Cut to the required size and shape. Secure with a cohesive or elastic support bandage such as NaturalintX Wrap. To avoid pressure points do not overlap or fold. At known pressure points a 'relief' hole should be made. Dressing should be changed frequently.*

CAUTION: *TO AVOID DANGER OF SUFFOCATION, KEEP PACKAGING OUT OF REACH OF BABIES AND CHILDREN. FOR ANIMAL USE ONLY. FOR EXTERNAL USE ONLY.*



Batch No: See Side Product code: 146
V2

For further information on using this
product contact: NAF, Wonastow Rd,
Ind Est West, Monmouth, NP25 5JA
call our Freephone Advice Line:
0800 373 106 or email info@naf-uk.com

