



V2

Batch No: See Side
Product code: 147

For further information on using
this product contact: NAF,
Wonastow Rd, Ind Est West,
Monmouth, NP25 5JA
call our Freephone Advice Line:
0800 373 106 or email info@naf-uk.com



ISO 9001
Quality Assured Firm



Wrap

10cm X 4.5m



NaturalintTMX



Wrap

10cm X 4.5m

DIRECTIONS FOR USE:

Apply the wrap at approximately 50% stretch.

Do not apply too tightly as this may impair circulation.

Apply last wrap with no tension.

Press firmly into place.

Wrap should be removed and re-applied daily.

Wrap should not be applied directly to open wounds.

CAUTION:

This product contains natural rubber latex which may cause allergic reactions in sensitive individuals.

**TO AVOID DANGER OF SUFFOCATION,
KEEP PACKAGING OUT OF REACH OF
BABIES AND CHILDREN. FOR ANIMAL
USE ONLY. FOR EXTERNAL USE ONLY.**